

No: 21195-202

BiClamp® 201 T



- Tissue structures can be coagulated over larger areas
- The individual dissection of vessels is usually not necessary
- Effective and rapid coagulation due to the optimal synchronization of the instrument and the VIO generator
- Highly economical: BiClamp is reusable; the procedure saves time and suture material
- AUTO STOP is user-friendly and offers optimal results
- When the instrument is plugged in, Plug & Play automatically adjusts the specific parameters

Dimensions and weight

Length	200 mm
Cable length	with connecting cable 4 m

Basic features

Shape / Forming	bent 18°
Connection	{MF-2-Stecker}

Surface

Branches	smooth
----------	--------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	0.25 kVp