

No: 10140-100

VIO[®] 300 D



- Automatic output dosage for all regulative technologies
- Precise power output adjusted to the individual indication
- New and improved CUT and COAG functions
- Variable configuration with socket modules
- Upgrades for specific medical specialties are possible at any time
- Simple, interactive and safe operation using the TFT display and FocusView
- Individually configurable programs in simple text: switch on the unit and begin working
- Plug and play: plug in the instrument, start work
- ReMode function: "remote control" via the handle or foot switch straight from the sterile field
- Assignment of the footswitch: the display makes it possible to assign the footswitch or even just one footswitch pedal to the output socket of choice

Power output

Maximum CUT output	300 watt at 500 ohm
Safety system	NESSY
Frequency	350 kHz
Maximum COAG output	up to 200 watt

Mains connection

Power frequency	50 / 60 Hz
Power input during stand-by	40 Watts
Power input during max RF output	500 watt / 920 VA
System voltage	100-120 / 220-240 Volt
Potential equalization connection	yes

Dimensions and weight

Width x Height x Depth	410 x 165 x 380 mm
Weight	9,5 kg

Additional product information

Packaging unit	1 Piece
Display	TFT-Display and FocusView

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF
Category acc. to UMDNS	11-490