


No: 20321-009

Evacuation tubing



Dimensions and weight

Length	2.7 meter
--------	-----------

Basic features

Compatibility	connection Ø 22 mm
---------------	--------------------

Connection	with connection Ø 22 mm
------------	-------------------------

Cleaning and Reprocessing

Sterilization in autoclave	max 134 °C
----------------------------	------------

Cleaning in washer disinfectant	max 95 °C
---------------------------------	-----------

Standards

Classification acc. to 93/42 EWG	Ila
----------------------------------	-----

Protection category acc. to EN 60 601-1	
---	--

Type acc. to EN 60 601-1	
--------------------------	--

Additional product information

Packaging unit	1 Piece
----------------	---------