

No: 20191-335



Coagulation electrode

Dimensions and weight

| | |
|----------------|--------|
| Shaft diameter | 2 mm |
| Length | 110 mm |

Basic features

| | |
|-----------------|----------------|
| Shape / Forming | tip angled 90° |
|-----------------|----------------|

Cleaning and Reprocessing

| | |
|---------------------------------|------------|
| Cleaning in washer disinfectant | max 95 °C |
| Sterilization in autoclave | max 134 °C |
| Single-use product | yes |

Standards

| | |
|---|-----|
| Classification acc. to 93/42 EWG | IIb |
| Protection category acc. to EN 60 601-1 | |
| Type acc. to EN 60 601-1 | |

Additional product information

| | |
|--------------------------|---------|
| Max. electrical capacity | 4 kVp |
| Packaging unit | 5 Piece |