

No: 20190-148

## Electrosurgical pencil for nerve test with 1 button, for silver-ball electrode



### Dimensions and weight

Cable length	with connecting cable 4 m
Shaft diameter	4 mm

### Basic features

Compatibility	NT 2 / ICC 350 ZMK
---------------	--------------------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Max. electrical capacity	0.25 kVp
Packaging unit	1 Piece