


No: 20191-004

## Sterilizing insert for 8 electrodes



- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines

### Dimensions and weight

Dimensions	66 x 17 x 46 mm
Shaft diameter	4 mm

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 134 °C

### Standards

Classification acc. to 93/42 EWG
Protection category acc. to EN 60 601-1
Type acc. to EN 60 601-1

### Additional product information

Packaging unit	1 Piece
----------------	---------