



No: 20195-014

Bipolar forceps Classic, angled

Dimensions and weight

| | |
|-----------------|-----------------|
| Length | 165 mm |
| Tip | tip 1 mm, blunt |
| Gripping length | 6 mm |

Basic features

| | |
|-----------------|--------|
| Shape / Forming | angled |
|-----------------|--------|

Cleaning and Reprocessing

| | |
|---------------------------------|------------|
| Cleaning in washer disinfectant | max 95 °C |
| Sterilization in autoclave | max 138 °C |

Standards

| | |
|---|-----|
| Classification acc. to 93/42 EWG | IIb |
| Protection category acc. to EN 60 601-1 | |
| Type acc. to EN 60 601-1 | |

Additional product information

| | |
|--------------------------|---------|
| Max. electrical capacity | 0.5 kVp |
| Packaging unit | 1 Piece |