

No: 20195-066



Bipolar forceps Classic, short angled

Dimensions and weight

Tip	tip 0.7 mm, fine
Length	165 mm
Gripping length	5 mm

Basic features

Shape / Forming	short-angled
-----------------	--------------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Max. electrical capacity	0.5 kVp
Packaging unit	1 Piece