

No: 20321-028

## Electrosurgical pencil with 2 buttons



### Dimensions and weight

Cable length with connecting cable 3 m

### Basic features

Compatibility VIO, ICC, ACC

Connection International

### Cleaning and Reprocessing

Sterile packaging E0

Single-use product yes

### Standards

Classification acc. to 93/42 EWG IIb

Protection category acc. to EN 60 601-1

Type acc. to EN 60 601-1

### Additional product information

Packaging unit 25 Piece

Max. electrical capacity 4.5 kVp