

No: 20195-148

## Insert for BiClamp<sup>®</sup> LAP forceps, fenestrated



### Dimensions and weight

Needle length	340 mm
Shaft diameter	5 mm
Length	340 mm

### Basic features

Shape / Forming	fenestrated
-----------------	-------------

### Surface

Branches	semi-deep
----------	-----------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	1 Piece
----------------	---------