

No: 21196-119

Bipolar connecting cable



Dimensions and weight

Cable length	{4,50 m}
--------------	----------

Basic features

Compatibility	for WOLF resectoscopes
Connection	{MF-U-Stecker}

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	I
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	1 kVp