

No: 20194-086

Return electrode cable



Dimensions and weight

Cable length	5 m
--------------	-----

Basic features

Compatibility	VIO, ICC, ACC, non-ERBE units
Connection	International

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	I
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	0.5 kVp