

No: 20190-166

Electrosurgical pencil with 2 buttons



Dimensions and weight

Shaft diameter	2.4 mm
Cable length	with connecting cable 3 m

Basic features

Compatibility	VIO, ICC, ACC, non-ERBE units
Connection	International

Cleaning and Reprocessing

Sterile packaging	E0
Single-use product	yes

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	25 Piece
Max. electrical capacity	5 kVp