

No: 20190-163

## Electrosurgical pencil with rocker switch



### Dimensions and weight

Shaft diameter	2.4 mm
Cable length	with connecting cable 3 m

### Basic features

Compatibility	VIO, ICC, ACC, non-ERBE units
Connection	International

### Cleaning and Reprocessing

Sterile packaging	E0
Single-use product	yes

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	25 Piece
Max. electrical capacity	5 kVp