

No: 21191-116

Tungsten needle electrode

- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines



Dimensions and weight

Shaft diameter	2.4 mm
Length	40 mm
Dimensions	ø 0,5 x 3 mm

Basic features

Shape / Forming	straight
Compatibility	VIO, ICC, ACC
Flexibility	rigid
Colour	black

Surface

Insulation	insulated
------------	-----------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	1 kVp