

No: 21191-115

Needle electrode



- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines
- The entire range available with both 4 mm as well as 2.4 mm shafts

Dimensions and weight

Shaft diameter	2.4 mm
Length	35 mm
Dimensions	ø 0,8 x 20 mm

Basic features

Shape / Forming	bent
Compatibility	VIO, ICC, ACC
Flexibility	rigid
Colour	black

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	5 Piece
Max. electrical capacity	4.3 kVp