

No: 21191-115

Needle electrode



- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines
- The entire range available with both 4 mm as well as 2.4 mm shafts

Dimensions and weight

| | |
|----------------|---------------|
| Shaft diameter | 2.4 mm |
| Length | 35 mm |
| Dimensions | ø 0,8 x 20 mm |

Basic features

| | |
|-----------------|---------------|
| Shape / Forming | bent |
| Compatibility | VIO, ICC, ACC |
| Flexibility | rigid |
| Colour | black |

Cleaning and Reprocessing

| | |
|---------------------------------|------------|
| Cleaning in washer disinfectant | max 95 °C |
| Sterilization in autoclave | max 138 °C |

Standards

| | |
|---|-----|
| Classification acc. to 93/42 EWG | IIb |
| Protection category acc. to EN 60 601-1 | |
| Type acc. to EN 60 601-1 | |

Additional product information

| | |
|--------------------------|---------|
| Packaging unit | 5 Piece |
| Max. electrical capacity | 4.3 kVp |