

No: 21191-107

## Spatula electrode, flexible

- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines



### Dimensions and weight

Shaft diameter	2.4 mm
Length	45 mm
Dimensions	2 x 6 mm

### Basic features

Shape / Forming	straight
Compatibility	VIO, ICC, ACC
Flexibility	flexible
Colour	black

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	5 Piece
Max. electrical capacity	4.3 kVp