


No: 20132-225

## FiAPC Probe 2200 C



- No preparation of the probe or cable required
- Lower costs per APC application
- No gas-filter replacement
- Plug and play: plug in the instrument, start work
- More safety through prevention of cross contamination

Dimensions and weight	
Dimensions	ø 2,3 mm
Length	2.2 meter
Basic features	
Flexibility	flexible
Beam forms	C = Circumferential Beam
Shape / Forming	C
Compatibility	{APC 300, APC 2, APC 3}
Cleaning and Reprocessing	
Single-use product	yes
Sterile packaging	E0
Standards	
Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	
Additional product information	
Max. electrical capacity	5 kVp
Packaging unit	10 Piece