

No: 20195-168

Bipolar forceps Classic, bayonet



Dimensions and weight

Tip	tip 1.2 mm, blunt
Gripping length	5 mm
Length	250 mm

Basic features

Shape / Forming	bayonet, angled downwards
-----------------	---------------------------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	0.5 kVp