

No: 21191-067



Tungsten snare electrode

- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines

Dimensions and weight

| | |
|----------------|----------|
| Dimensions | 5 x 5 mm |
| Shaft diameter | 4 mm |
| Length | 125 mm |

Basic features

| | |
|-----------------|---------------|
| Shape / Forming | straight |
| Compatibility | VIO, ICC, ACC |
| Flexibility | rigid |
| Colour | black |

Cleaning and Reprocessing

| | |
|---------------------------------|------------|
| Cleaning in washer disinfectant | max 95 °C |
| Sterilization in autoclave | max 138 °C |

Standards

| | |
|---|-----|
| Classification acc. to 93/42 EWG | IIb |
| Protection category acc. to EN 60 601-1 | |
| Type acc. to EN 60 601-1 | |

Additional product information

| | |
|--------------------------|---------|
| Packaging unit | 1 Piece |
| Max. electrical capacity | 1.5 kVp |