

No: 20150-060

## HybridKnife® T-Type I-Jet



- No change of instrument necessary, thus shorter operating times
- Separation medium can be replenished quickly at any time
- Mechanical and thermal protective function of water cushion is preserved throughout the entire process of resection
- Good view of the operating site: risk of bleeding is minimized (blood vessels are compressed by the water cushion)
- Waterjet pressure can be adjusted depending on lesion, tissue layer and area of application

### Dimensions and weight

Dimensions	ø 2,3 mm
Tip	T-Type
Needle length	5 mm
Length	1.9 meter

### Basic features

Compatibility	VIO, ERBEJET 2
Flexibility	flexible
Connection	International (3-Pin)

### Cleaning and Reprocessing

Sterile packaging	E0
Single-use product	yes

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Max. electrical capacity	2.5 kVp
Packaging unit	5 Piece