

No: 21191-459

## Coated Spatula Electrode



### Dimensions and weight

Dimensions	2,3 x 19 mm
Shaft diameter	2.4 mm
Length	50 mm

### Basic features

Shape / Forming	straight
Colour	blue

### Surface

Coating	{PTFE-Beschichtung}
---------	---------------------

### Cleaning and Reprocessing

Sterile packaging	EO
Single-use product	yes

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	10 Piece
Max. electrical capacity	4 kVp