

No: 21191-181

Large surgical set, ø 2.4 mm



Dimensions and weight

Shaft diameter	2.4 mm
----------------	--------

Basic features

Compatibility	VIO, ICC, ACC
Flexibility	rigid
Colour	black

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	4.3 kVp