

No: 10402-000

## ERBECRYO<sup>®</sup> 2



- Convenient operation with plug and play: ERBECRYO 2 automatically sets the appropriate parameters for every probe
- Even more convenience: all important settings appear on the display: - Effect levels - Timer function
- Activation using the footswitch
- Reproducible freeze performance from the first time the probe is used through to the last

Mains connection	
System voltage	100-120 / 220-240 Volt
Power frequency	50 / 60 Hz
Potential equalization connection	yes
Dimensions and weight	
Width x Height x Depth	410 X 130 X 370 mm
Weight	6,7 kg
Additional product information	
Packaging unit	1 Piece
Programmspeicherplätze	{10}
Standards	
Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	I
Type acc. to EN 60 601-1	CF